

BLANCO

Portable suction pump 12V 16 l/min with battery



BLANCO PORTABLE ASPIRATOR 12V 16 l/min with battery is a device working with 230V ~ AC 12V DC mixed power supply.

Easily trasportable from one hospital ward to another, particularly suitable for tracheotomized patients, minor surgical applications and post – operative therapy at home to be used for the nasal aspiration, oral aspiration, tracheal aspiration of the body liquids (mucus or catarrh) in the adult or in the children.

Easily portable equipment designed for virtually continuous use, thanks to an electronic system for managing the electricity supply with in addition to monitoring the functioning of the internal battery, also enables the equipment to be used at a voltage of 230 V ~ even when the battery is being re-charged.

A series of LED on the front control panel signals whether the machine is operating, as well as the charge and re-charge status of the batteries.

Made of highly heat - resistant, electrically insulated plastic material in conformity with the latest European safety standards.

Supplied with a polycarbonate autocavable jar complete with overflow valve.

Equipped with aspiration regulator, located on the front panel.

Technical details

Classification UNI EN ISO 10079-1	HIGH VACUUM / LOW FLOW
Security class Feeding	Class II type B 230V – 50 Hz – 12Vdc
Power consumption	50 VA
Fuse	F1 x 630 mA 250V
Max aspiration (without jar)	-75 kPa (-0.75 bar)
Minimum aspiration	Less -40 kPa (-0.40 bar)
Max flow aspiration (without jar)	16 L/min
Weight	4.9 kg
Dimensions	370x220x210 mm
Battery holding time	45/50 minutes
Battery time charge	120/150 minutes

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
E-mail: info@spencer.it Web: <http://www.spencer.it>

This document has been automatically generated by www.spencer.it on 2010/01/7 12:12:19, and is referred to the last production lot. The informations furnished have been last updated on 2009/12/01. To verify the correctness of this document, please contact our Quality Dept. at quality@spencer.it

SC10002B

Blanco portable
suction device

